

# Quantification of mRNA Expression using DNA-based Standard Curve RT-qPCR Methods

Milena Blaga<sup>1\*</sup>, Neil Henderson<sup>2</sup>, Wayne Petrie<sup>1</sup>, Kirsty Edwards<sup>1</sup>, Gilta Jaeckel<sup>1</sup>, Alexander Gunn<sup>1</sup>, Stephen Cranwell<sup>1</sup>, Alice Ciobanu<sup>1</sup>, Louise Baillie<sup>1</sup>, Craig Stovold<sup>3</sup>, Keith Sutton<sup>1</sup>, Stephen Williams<sup>1</sup>, Richard Heideman<sup>1</sup>

<sup>1</sup> Department of Immunobiology, Charles River Laboratories Edinburgh Ltd., UK.  
<sup>2</sup> Translational Biomarkers and Bioanalysis, Clinical Pharmacology and Safety Sciences, R&D, AstraZeneca, Gothenburg, Sweden.  
<sup>3</sup> Translational Biomarkers and Bioanalysis, Clinical Pharmacology and Safety Sciences, R&D, AstraZeneca, Cambridge, UK.  
\* E-mail: [milena.blaga@crl.com](mailto:milena.blaga@crl.com)

## 1 INTRODUCTION

Gene expression analyses are bioanalytically relevant for quantifying expression strength of vector-based gene therapeutics (TK/PK), as well as for establishing efficacy (PD) of oligonucleotide-based therapeutics (Figure 1). Reverse transcription quantitative PCR (RT-qPCR) is often the method of choice for performing gene expression assessments, with alternative methodologies including branched DNA assays (QuantiGene), NanoString or droplet digital PCR (ddPCR). RT-qPCR methods vary in technical considerations and analytical design:

- one-step or two-step
- oligo dT, random, or gene specific primers
- SYBR or probe-based
- singleplex or multiplex reactions
- relative to a DNA/RNA calibration standard
- relative to reference genes ( $\Delta\Delta C_t$ , Livak and Schmittgen, 2001)

We have developed a simple and efficient workflow for validating RT-qPCR assays for mRNA quantification using a single-step, DNA standard curve-based RT-qPCR method.




Figure 1: Gene expression modulation in disease and in response to therapeutic intervention.

## 2 MATERIALS AND METHODS

RNA Extraction performed using a phenol-chloroform based method (TRIzol Reagent, Invitrogen)



Figure 2a: Workflow for the preparation of purified RNA from tissues.

DNase-treatment using the DNA-free DNA Removal Kit (Ambion)  
RNA quantification using fluorescence-based Quant-iT RiboGreen RNA Assay Kit (Invitrogen)




Figure 2b: Workflow for RNA processing prior to RT-qPCR analysis.

RT-qPCR in a single-step reaction using QuantiTect Multiplex RT-PCR Kit (Qiagen). Detection of target can be performed using custom-designed primer-probe sets or off-the-shelf assays. For this case-study, a Bio-Rad assay and DNA template ( $2 \times 10^7$  copies/ $\mu$ L) were used for Gene X quantification in mouse liver and placenta. Data was captured on a QuantStudio 7 Flex System and analysed in Watson LIMS.

## 4 RESULTS




Figure 4a: Calibration standards linearity. Lot to lot variation impact on precision (Ct signal) but not on accuracy (concentration).

Assay Performance

	Accuracy (RE)	Precision (CV)
Intra-Assay	$\pm 1.4\%$	$\leq 0.6\%$
Inter-Assay	$\pm 0.7\%$	$\leq 0.6\%$

Table 1: Intra-/Inter-assay accuracy and precision results based on single reference standard lot numbers.

Gender Differences in Liver Gene X Levels




Figure 4b: Non-treated liver gene expression in 25 male and 25 female mice. Assessment of control levels and gender differences.

Gene X Levels in Placenta






Figure 4c: Non-treated placenta gene expression in 10 female mice at 2 concentrations, showing linear fold-change. Reduction of input RNA simulates mRNA knockdown expected in treated animals.

## 3 METHOD VALIDATION


Reference Standard and RNA Linearity



Limit of quantification/detection (LOQ/LOD)



Inter-/Intra-Assay Accuracy and Precision



RNA Levels in Individual Control Samples




Figure 3: Plate map templates for characterising assay performance. Cal – Calibration Standard. WB – Water Blank. No RT – No Reverse Transcriptase. QC – Quality Control (Upper Limit, High, Mid, Low, Lower Limit). S – Sample.

## 5 CONCLUSION

- Application criteria:
- Linearity is assessed in 2 samples for each matrix. LOQ/LOD is assessed once with DNA standard.
  - Intra-/inter-assay accuracy and precision is assessed on 6 occasions by at least 2 analysts.
  - Control sample analysis should be performed for each matrix and gender, at 2 RNA concentrations.
  - Calibration standards should be analysed in triplicate wells; QCs, samples and blanks in duplicate wells.
  - At least 6 calibrators, 66% of QCs and one QC at each level should meet the acceptance criteria.
  - The coefficient of determination  $r^2$  should be  $\geq 0.990$ . The amplification efficiency should be 90-110%.
  - Precision (%CV) of Ct values should be  $\leq 3\%$ . Accuracy (%RE) of log copies should be within  $\pm 10\%$ .
  - Water blanks and No RT controls containing RNA should be  $< LOD$ .

A single-step RT-qPCR is advantageous over a two-step protocol based on: reduced potential for contamination, lower variability due to single-reaction setup, lower effects from RNA integrity due to gene-specific primer-probe sets (Fleige and Pfaffl, 2006), as well as time-efficient and simple setup, reducing the likelihood of process errors. The minimum-required testing parameters described in this workflow allow gene expression assays to be validated in a straightforward and reproducible manner, using any target sequence and custom primer-probe sets.

References:  
Livak and Schmittgen, Methods, 2001,25(4):402-8, DOI: 10.1006/meth.2001.1262  
Fleige and Pfaffl, Molecular Aspects of Medicine, 2006, 27:126-139, DOI: 10.1016/j.mam.2005.12.003